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Contents

Overview	5
Safety instructions	7
Unpacking / scope of delivery	9
Installation	11
Operation	13
Switching the device on/off	13
Display	15
Changing the temperature Setpoint	17
Changing the measurement unit	19
What to do, if...	21
Maintenance and disposal	23
Cleaning	23
Disposal	23
Calibration service	23
Technical data	25
Approvals	27

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
Overview

The devices GFX 460, GFX 460 G and GFX 460-B are waterproof precision regulating thermometers for temperatures between -50 °C to +300 °C (-58 °F to +572 °F). They are capable of controlling the heating of laboratory devices such as magnetic stirrers, heating stirrers etc.

Overshooting of the temperature is minimized by a fuzzy control. By means of this type of control the heating is switched off before reaching the desired value.

The GFX 460 device is fitted with a stainless steel probe (Ø 3 mm) and the GFX 460 G device with a glass probe (Ø 7 mm) for use with especially aggressive substances. Different Pt100 probes can be plugged on GFX 460-B.

Device elements:

1. Display field (LCD) with background lighting
2. Heating display (LED green/red)
3. Button ▼ (reduce desired temperature)
4. SET button
5. Temperature probe
6. connection to stirrer/current supply
7. Button ▲ (increase desired temperature)

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Safety instructions

Do not use the device in potentially explosive areas. Mortal danger!

Do not place the device in an environment that has temperatures exceeding 70 °C. Potential danger of explosion of the rechargeable lithium battery!

Do not subject the device to microwave radiation. Potential danger of explosion of the rechargeable lithium battery!

Do not touch the probe with your bare hand immediately after measuring a high temperature. Risk of burns!

In order to avoid damage to the device, protect it from direct contact with water. Observe the protection class of the device (IP65).

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Unpacking / scope of delivery

Check that the contents are in perfect condition after receiving the device package. Also check that the contents of the device package comply with your order.

If you have any reason for complaint, do not hesitate to contact us. You will find our contact data on the back of this manual.

Scope of delivery

The delivery includes these parts:

- This user manual
- GFX 460 Measuring Device with stainless steel probe, GFX 460-G with glass probe or GFX 460-B with pluggable Pt100 probes
- a calibration service card (see page 23)
- device fixture

Possible accessory items

In addition various accessories may be included in the deliver:

Description	Identifier
Mains adaptor	AX 410
Relay box	AX 400

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
Installation

The device needs a stable voltage supply of 9...12 V. This is provided for many electrical laboratory devices by means of a 5 pin DIN plug socket.

For this reason the device has a connector cable with a DIN diode plug that also contains the heating control cable. This connector cable is designed for direct connection to the heating plates from the various leading manufacturers.

If needed you can find the plug configuration in the adjacent.

You will need the relay box RB 1691 that is available as an accessory in order to connect heating plates or heating mantles that do not have a DIN plug socket.

- Position the device at the desired point.
- Put the DIN diode plug on the connector cable into the socket of the laboratory device to be used.

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
Operation

Switching the device on/off

Current for the measuring device is supplied from the connected stirring device (optionally: by using the mains adapter provided). The measuring device starts up and is ready for use in seconds (self-testing) as soon as the stirring device is switched on or the mains adapter is plugged in.

The measuring device is switched off simultaneously as soon as the stirring device is switched off (or mains adapter is plugged out).

The display background lighting and the heating display (LED, see page 5) continue to light up for as long as the measuring device is supplied with current.

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

Display

The display field (LCD) comprises the following display elements:

- 1 Tendency indicator ▲ / ▼
- 2 Current measured value
- 3 Desired value/text
- 4 Measurement unit (°C or °F)

The measured value presentation

During operation the device displays the temperature continuously in either °C or °F depending on the setting (see page 19).

The device updates the value displayed in one second cycles. When the temperature is rising or falling a tendency display appears to the left of the current temperature value.

The heating display

When the controls are set to heat, the device switches the heating contact by means of which a laboratory device is controlled to “on”. The color of the light emitting diode under the display field now changes from green to red. Simultaneously “HEAT” appears in the display field below the current measured value.

The device switches the heating contact to “off” again shortly before the measured value reaches the desired value. The light emitting diode changes back from red to green. The “HEAT” display fades out.

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

Changing the temperature Setpoint

- Press the SET button to change the temperature setpoint.

The temperature setpoint valid until now appears below the current measured value.


- Use the ▲ and ▼ buttons to change the temperature setpoint.
 - One short press changes the setpoint by 1°.
 - For setpoint changes of a greater range press and hold the button for a longer period.
- Press the SET button again to save the temperature setpoint.

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Changing the measurement unit

You can set whether the device should measure the temperature in °Celsius or °Fahrenheit.

Call up the User Menu to make this setting. Now proceed as follows:

- Prerequisite: the device is switched off (switch of stirring device or disconnect cable if necessary).
- Press and hold the SET button and then switch on current (switch on stirring device or connect cable or plug in optional mains adapter).
- Press and hold button until "USER" appears on display.

After a short wait "UNIT" appears on display.

- To change the measurement unit from °Celsius or °Fahrenheit or back again press SET.

The change is saved at once.

- To leave the User Menu switch off the device by switching off the connected stirring device, disconnecting the cable or pulling the mains adapter plug.

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Maintenance and disposal

Cleaning

Clean the device with a slightly dampened cloth.

Never use solvents (such as acetone) for cleaning because these may attack the plastic.

Disposal

If the device becomes no longer fit for purpose, it must be disposed of in a suitable, environmentally-friendly manner.

Never dispose of the device in the domestic garbage.

Calibration service

In order to guarantee a high measuring precision, the thermometer must be calibrated every year. Therefore we offer you a calibration service.

- Fill in the enclosed service card.
- After one year we will recall your thermometer for calibration.

Send us your devices and we will send them after calibration back.

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Technical data

Measurement	Temperature (°C/°F)
Probe	
GFX 460	Pt1000
GFX 460-G	Pt1000
GFX 460-B	Pt100, 4-wire technologie
Measuring range	-50 °C...+300 °C -58 °F...+572 °F
Measuring precision	±1 °C (±1 °F)
Measurement interval	1s
Resolution of measuring signal	1 °C, 1 °F
Thermal time constant T99 flowing water	
GFX 460 (stainless steel probe) ca. 8 s	
GFX 460 G (glass probe) ca. 8 s	
Thermal time constant T99 still Wasser	
GFX 460 (stainless steel probe) ca. 10 s	
GFX 460 G (glass probe) ca. 90 s	
Operating temperature	-5 °C...+60 °C (+23 °F...+140 °F)
Storage	-30 °C...+70 °C (-22 °F...+158 °F)
Control type	fuzzy control
Protection type	IP65
Power supply external	9...12 V DC, 20 mA
Safety features	
Probe breakage recognition	
probe terminal recognition	
Safety circuit and control circuit to DIN 12878 classes 1 and 2	
Other characteristics	
Display background lighting	
LED indicator for control status	
Connection	5 pin diode plug DIN 45322 see illustration p. 10 for pin assignment
Dimensions	110 x 54 x 22 mm (L x W x H) without probe
Housing material	Acid resistant polyamide
Weight	approx. 150 g

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Approvals

According to the certificate of conformity, this product complies with the EMC directive.

Information, operating instructions and declarations of conformity can be found at www.ebro.com

Operation

Operation is subject to the following two conditions:

- (1) This device may not cause harmful interference, and
- (2) This device must accept any interference received, including interference that may cause undesired operation.

Changes or modifications of the device could void the user's authority to operate the equipment.

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Service-Adresse / Service Address / Adress du Service

Xylem Analytics Germany Sales GmbH & Co. KG
ebro
Am Achalaich 11
82362 Weilheim
Germany

Phone: +49.(0)841.954.78.0
Fax: +49.(0)841.954.78.80
Internet: www.ebro.com
E-Mail: ebro@xylem.com

.....

Hersteller / Producer / Fabricant

Xylem Analytics Germany GmbH
Am Achalaich 11
82362 Weilheim
Germany