

DAC Universal MK III and MK IV

SEMINARS AND VALIDATION



Seminars: Process validation with devices DAC Universal MK III and MK IV

We train you in our KompetenzCentrum in the application with the validation equipment directly at the DAC Universal.

Trainings in German language

May 8th - 9th 2019 and June 4th - 5th 2019
Ingolstadt, Germany

Validierung DAC Universal MK III + MK IV

Contact: +49 841 95478-0 or
<https://www.ebro.com/de/unternehmen/kompetenzcentrum/>

Equipment for validation DAC Universal MKIII and MKIV



SL 3301 Data logger set for thermal tests in the **DAC Universal**

(According to EN 13060 and EN ISO 15883)

Part No. 1250-3301



AL 3300 Thermal test on the **DAC MK III** and **DAC MK IV**

(According to EN 13060 and EN ISO 15883)

Part No. 1340-6052



AL 3301 - AL 3303 Testing of cleaning performance on **DAC MK III** and **DAC MK IV**

(According to EN ISO 15883 with test soil supplied by SMP)

AL 3301 **Part No. 1248-3301**

AL 3302 **Part No. 1248-3302**

AL 3303 **Part No. 1248-3303**



AL 3304 Connecting column, additionally required for DAC MK IV for thermal tests on the **DAC MK IV**

(According to EN ISO15883)

Part No. 1248-3304

Headquarters:

Xylem Analytics Germany Sales GmbH & Co. KG, ebro
Peringerstr. 10 · 85055 Ingolstadt, Germany
Phone +49 841 95478-0 · Fax +49 841 95478-80
Internet: www.ebro.com · E-Mail: ebro@xyleminc.com

Representatives worldwide:

www.ebro.com