

Temperature Data Logger EBI 12

for monitoring, process and routine control as well as temperature mapping

BLUE LINE

The EBI 12 data logger is distinguished by its robustness and versatility.

This logger was specially developed for temperature measurement in the range around -80 °C.

- High-quality stainless-steel housing
- Operating range from -90 °C to +60 °C
- High temperature accuracy of up to ±0.3 °C
- sterilizable

EBI 12-T102



Technical Data:

of

Standards and Directives:

- 2014/53/EU
 Radio Equipment Directive RED
- (EC) 1935/2004
 Materials and articles which come into contact with food

Туре	Description	Part No.
FBI 12-T102	Data Logger with internal sensor	1340-6636

