

Service order

For cold chain data loggers of the product series:
EBI 20, EBI 25, EBI 300, EBI 310, EBI 40

For internal use:

Document version: CCL.2022.02.EN



SI Analytics®



are brands of Xylem Inc. or one of its subsidiaries.

Your Order number:

Billing Address:

Your customer number (if known):
VAT number:
Company:
Address1:
Address2:
Address3:
ZIP-Code:
City:
Country:

Delivery Address (if different):

Number (if known):
VAT Number:
Company:
Address1:
Address2:
Address3:
ZIP-Code:
City:
Country:

Contact Information:

Name:
Phone:
eMail:

Certificate Address (if different):

Company:
Address1:
Address2:
Address3:
ZIP-Code:
City:
Country:

We herewith order the service selected on the following pages as well as the professional disposal of the replaced parts. The order is placed under the general terms and conditions of Xylem Analytics Germany Sales GmbH & Co. KG. <https://www.xylemanalytics.com/en/terms-and-conditions>. The prices do not include value-added tax or the costs of packaging, [...], transportation and insurance, if applicable. Delivery terms are valid according to INCOTERMS 2010.

Date

Name in letters

Signature

You can use this form to order the same service for one or more data loggers. If you like to order a different service for other loggers, please use a separate form.

Please pack them separately to possible other service orders to avoid mistakes.

This document is valid until December 31st 2022, version CCL.2022.02.EN



Service order

For Data Loggers of the product series EBI 20, EBI 25, EBI 300, EBI 310, EBI 40.

Your Order Number:

(Only necessary if you order more than one different services.)

Serial Number (s):

Model or product description:

Issue description:

(if any)

ISO Calibration:

ISO Calibration at two standard points for 78.00 € *¹)

ISO Calibration at three standard points for 112.00 € *¹)

ISO Calibration at four standard points for 145.00 € *¹)

ISO Calibration at five standard points for 179.00 € *¹)

ISO Calibration at two standard temperature points and two humidity points for 102.00 € *¹)

Acc. calibration acc. to DIN EN ISO/IEC 17025 (DAkKS) at calibration points of your choice:

Two temperature points for 195.00 €

Three temperature points for 242.00 €

Four temperature points for 283.00 €

Five temperature points for 324.00 €

Two temperature points and two humidity points for 469.00 €

Three temperature points and three humidity points for 510.00 €

Calibration points of your choice:

Warranty repair. (to verify the claims, please attach the appropriate copy of your purchase receipt)

Please provide cost estimate. Note: for fastest service, please select one of the standard service options

Select Desired Options (available for an additional fee):

As Found as Left calibration

Battery replacement

Probe replacement

*¹) please take a look at additional information on following pages for further information.

www.wtw.com

www.si-analytics.com

www.ebro.com

www.meinsberg.de

Xylem Analytics Germany Sales GmbH & Co. KG
Am Achalaich 11
82362 Weilheim, Germany
Amtsgericht München HRA 105345
USt-IdNr./VAT No.: DE815629128
Steuernummer: 168/183/00305

Bank: Deutsche Bank AG
BLZ: 700 700 10
Kto.-Nr.: 992420000
IBAN: DE11700700100992420000
BIC: DEUTDEMMXXX
WEEE Nr./No.: DE 29884659

Persönlich haftende Gesellschafterin: Xylem Analytics Germany GmbH
Sitz: Weilheim i.OB,
Registergericht: Amtsgericht München, HRB 145685
Geschäftsführer/Managing Directors:
Dr. Robert Reining · Ulrich Schwab ·
USt-IdNr./VAT No.: DE275429984

You can use this form to order the same service for one or more data loggers. If you like to order a different service for other loggers, please use a separate form.

Please pack them separately to possible other service orders to avoid mistakes.

This document is valid until December 31st 2022, version CCL.2022.02.EN



Additional Information

What does the service include?

ISO calibration: The device will be calibrated according to internal work orders and defined calibration methods. The measured values were compared to the etalon who is traceable to national standard.

Accredited calibration according to DIN EN ISO / IEC 17025: Calibration according the DIN EN ISO / IEC 17025 (DAkKS) requirements.

(1) What are our standard calibration points?

Product series	Logger types	Standard calibration points
EBI 20	EBI 20-T1, EBI 20-TE1	-20 °C, 0 °C
	EBI 20-TF	+20 °C, +60 °C
	EBI 20-TH1	0 °C, +20 °C 32.8 %rH, 75.4 %rH at 25 °C
EBI 25	EBI 25-T, EBI 25-TE	-20 °C, 0 °C
	EBI 25-TX with one probe	-80 °C, 0 °C, +134 °C
	EBI 25-TH	0 °C, +20 °C 32.8 %rH, 75.4 %rH at 25 °C
EBI 300	EBI 300, EBI 300 TE	-20 °C, 0 °C
	EBI 300 TH	0 °C, +20 °C 32.8 %rH, 75.4 %rH at 25 °C
EBI 310	EBI 310	-20 °C, 0 °C, +60 °C
	EBI 310 TE completet	Logger: -20 °C, 0 °C, +60 °C Probe: -80 °C, 0 °C, +134 °C
	TPX 220 probe	-80 °C, 0 °C, +134 °C
	EBI 310 TX without probe	-200 °C, -20 °C, 0 °C, +60 °C, +400 °C
	TPX 310 Messadapter without probe	-200 °C, 0 °C, +400 °C
	EBI 310 DI	-80 °C, -20 °C, 0 °C, +60 °C
	TPX 250 probe	-80 °C, 0 °C
	EBI 310 TH	0 °C, +20 °C 32.8 %rH, 75.4 %rH at 25 °C
EBI 40	EBI 40 without probe	-100 °C, 0 °C, +1,000 °C

You can use this form to order the same service for one or more data loggers. If you like to order a different service for other loggers, please use a separate form.

Please pack them separately to possible other service orders to avoid mistakes.

This document is valid until December 31st 2022, version CCL.2022.02.EN



What are the prices for ISO calibration?

Item Number	Description	Price
1030-2302	ISO Calibration at 2 temperature points	78.00 €
1030-2303	ISO Calibration at 3 temperature points	112.00 €
1030-2304	ISO Calibration at 4 temperature points	145.00 €
1030-2305	ISO Calibration at 5 temperature points	179.00 €
1030-2322	ISO Calibration at 2 temperature and 2 humidity points	102.00 €

What are the prices for the accredited calibration according to DIN EN ISO / IEC 17025 (DAkkS) calibration?

Item Number	Description	Price
1030-3302	Accredited calibration at 2 temperature points	195.00 €
1030-3303	Accredited calibration at 3 temperature points	242.00 €
1030-3304	Accredited calibration at 4 temperature points	283.00 €
1030-3305	Accredited calibration at 5 temperature points	324.00 €
1030-3306	Accredited calibration at 6 temperature points	371.00 €
1030-3322	Accredited calibration at 2 temperature and 2 humidity points	469.00 €
1030-3323	Accredited calibration at 3 temperature and 3 humidity points	510.00 €
1030-3324	Accredited calibration at 4 temperature and 3 humidity points	556.00 €
1030-3325	Accredited calibration at 5 temperature and 3 humidity points	597.00 €
1030-3326	Accredited calibration at 6 temperature and 3 humidity points	639.00 €