

Software

ebro® offers exactly the software you need:

- Evaluation software for any applications:
Winlog.pro
- Evaluation software for EBI 25 data loggers:
Winlog.web
- Evaluation software for pharmaceutical and medical applications:
Winlog.med and Winlog.validation

Software/Features	Winlog.pro	Winlog.med	Winlog.validation	Winlog.web
Event Triggered Recording	●	●	●	●
Script-Calculations	●			●
System-Scripts	●			
Picture Manager	●		●	
Measure in Charts	●			
Cursor	●	●	●	
Realtime Calculations	●	●	●	●
Ranges	●	●	●	
Range-based Calculations	●	●	●	
Statistics per Range	●	●	●	
Relative Time Axis	●			
Configurations	●	●(2)	●(2)	
Firmware-Update	●	●	●	●
Import	●(1)			
Calibration	●	●	●	
Automatic File Name Generation	●			
21 CFR Part 11	●	●	●	●
User Administration	●	●	●	●
Audit-Trail	●	●	●	●
Advanced Chart Features	●			
Multi Document Support	●	●	●	
Export (Excel, PDF)	●	●	●	●
Customizable Company Logo	●	●	●	●
Wireless Support	●	●	●	●
2D Placement	●		●	●
3D Placement			●	
Routine Check	●(3)	●	●	
Validation	●(3)		●	
Unit Administration	●	●	●	●
Split Measurements	●	●	●	
Advanced Alarm-Management (Zones)	●(5)			
IQ/OQ	●		●	●

(1) Winlog.validation
(2) Template based

(3) Manual
(4) Since V2.5
(5) EBI 310 only

(6) Since V2.6
(7) Since V3.3
(8) Since V3.5

(9) Since V2.8
(10) Since V3.7



Supported Logger Types	Winlog.pro	Winlog.med	Winlog.validation	Winlog.web
EBI 10	●	●	●	
EBI 11	●	●	●	
EBI 12	●(9)	●(10)	●(10)	
EBI 16		●(8)	●(8)	
EBI 20	●	●(7)	●(7)	
EBI 25				●
EBI 40	●	●(7)	●(7)	
EBI 100	●	●	●	
EBI 310	●(6)	●(7)	●(7)	
CT 830	●	●	●	
PHT 830	●	●	●	

System Requirements	Winlog.pro	Winlog.med	Winlog.validation	Winlog.web
Windows 8 / 32 bit	●	●	●	●
Windows 8 / 64 bit	●	●	●	●
Windows 10 / 32 bit	●	●	●	●
Windows 10 / 64 bit	●	●	●	●
Windows 11 / 64 bit			●	●
Memory	≥ 1GB	≥ 2 GB (32bit) / ≥ 4 GB (64bit)	≥ 2 GB (32bit) / ≥ 4 GB (64bit)	≥ 4 GB
Hard Disc Memory	≥ 1GB	≥ 4 GB	≥ 4 GB	≥ 20 GB
DVD Drive	●	●	●	●
Screen Resolution	≥ 1280x768	≥ 1280x768	≥ 1280x768	≥ 1280x768
Processor	Dual Core 1.6 GHz+	Dual Core 1.6 GHz+	Dual Core 1.6 GHz+	Dual Core 2 GHz+

Market Overview	Winlog.pro	Winlog.med	Winlog.validation	Winlog.web
Food	●			●
Industry	●		●	●
Pharmaceutical	●		●	●
Medical	●	●	●	●

Available Languages	Winlog.pro	Winlog.med	Winlog.validation	Winlog.web
English	●	●	●	●
French	●	●	●	●
Italian	●	●	●	●
Spanish	●	●	●	●
Chinese	●	●	●	●
Japanese	●	●	●	
Korean	●	●	●	
Portuguese		●	●	
Greek		●	●	
Czech	●	●	●	●
Swedish	●	●	●	
Dutch	●	●	●	
German	●	●	●	●
Turkish		●	●	
Polish		●	●	
Norwegian		●	●	

Software

ebro® offers exactly the software you need:

- Evaluation software for any applications:
Winlog.basic and Winlog.pro
- Evaluation software for EBI 25 data loggers:
Winlog.web and Winlog.wave

Installation Qualification and Operation Qualification documents (IQ/OQ) are available for the Winlog.web and Winlog.pro softwares. They are written in English and saved in an editable file format (MS Word and/or MS Excel). They contain instructions on how to check the functionality of the software features. The instructions can be expanded on demand. Our support personnel is happy to help you executing the qualification.

Software/Features	Winlog.basic	Winlog.pro	Winlog.wave	Winlog.web
Event Triggered Recording		●		●
Script-Calculations		●		●
System-Scripts		●		
Picture Manager		●		
Measure in Charts		●		
Cursor		●		
Realtime Calculations		●		●
Ranges		●		
Range-based Calculations		●		
Statistics per Range		●		
Relative Time Axis		●		
Configurations		●		
Firmware-Update		●		●
Import		●(1)		
Calibration		●		
Automatic File Name Generation		●		
21 CFR Part 11	●	●	●	●
User Administration	●	●	●	●
Audit-Trail	●	●	●	●
Advanced Chart Features		●		
Multi Document Support		●		
Export (Excel, PDF)	●	●	●	●
Customizable Company Logo	●	●	●	●
Wireless Support		●	●	●
2D Placement		●	●	●
Unit Administration	●	●		●
Split Measurements		●		
Advanced Alarm-Management (Zones)	●(2)	●(2)		
IQ/OQ		●		●

(1) From Winlog.basic

(2) EBI 310 only

(3) Since V2.5

(4) Since V2.6

(5) Since V2.63



Supported Logger Types	Winlog.basic	Winlog.pro	Winlog.wave	Winlog.web
EBI 20	●	●		
EBI 25			●	●
EBI 40		●		
EBI 300	●(3)	●(3)		
EBI 310	●(4)	●(4)		
EBI 330	●(5)	●(5)		
CT 830		●		
PHT 830		●		

System Requirements	Winlog.basic	Winlog.pro	Winlog.wave	Winlog.web
Windows 8 / 32 bit	●	●	●	●
Windows 8 / 64 bit	●	●	●	●
Windows 10 / 32 bit	●	●		●
Windows 10 / 64 bit	●	●		●
Memory	≥ 1 GB	≥ 1 GB	≥ 1 GB	≥ 4 GB
Hard Disc Memory	≥ 1 GB	≥ 1 GB	≥ 1 GB	≥ 20 GB
CD/DVD Drive	●	●	●	●
Screen Resolution	≥ 1280x768	≥ 1280x768	≥ 1280x768	≥ 1280x768
Processor	Dual Core 1,6 GHz+	Dual Core 1,6 GHz+	Dual Core 2 GHz+	Dual Core 2 GHz+

Market Overview	Winlog.basic	Winlog.pro	Winlog.wave	Winlog.web
Food	●	●	●	●
Industry	●	●	●	●
Pharmaceutical		●		●

Available Languages	Winlog.basic	Winlog.pro	Winlog.wave	Winlog.web
English	●	●	●	●
French	●	●	●	●
Italian	●	●	●	●
Spanish	●	●	●	●
Chinese	●	●	●	●
Japanese	●	●		
Korean		●		
Czech	●	●		●
Swedish	●	●		
Dutch		●		
German	●	●	●	●